



## Stomacher® 4500 BAM



The robust Stomacher® 4500 is ideal for the BAM large volume testing procedures.

Expanding our range of Blenders focused on food laboratories requirements, the Stomacher® 4500 BAM is the perfect tool for extraction from 375g composite samples.

The Stomacher technique has become the world standard in Sample preparation. The heavy duty Stomacher 4500 BAM complements the Stomacher 400 Circulator

and provides the food laboratories with a high performance, high quality blender capable of homogenizing samples up to 4500ml.

### Why the Stomacher® 4500 BAM?




Ideal for 375g samples

Efficient Large Volume Stomacher

Volume sample size:  
1000 – 4500ml

Programmable Settings

3 year Warranty

Machine Model	Paddle Type	Sample weight*	Diluent Volume*
 <b>Stomacher 4500</b>		<b>375g</b>	<b>3375ml</b>
Bags/Accessories	Standard Bags	Bag Rack	Bag Clip
	BA6042 BA6042/HD BA6042/CLR	BA6094	BA6097

For further information on Seward products please contact your local Seward distributor.

### Seward (UK)

Europe, Middle East, Indian Subcontinent, Africa, Asia, Pacific, ROW

**+44 (0)1903 524600**

info@seward.co.uk

### Seward (USA)

North America, South America

**+1 631 337 1808**

info@sewardusa.com

\* as per the FDA BAM protocol protocol listed 20 September 2014 : [www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods](http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods)