



Stomacher® 80 Biomaster



Small tissue processing solutions for clinical and life science microbiology

The Stomacher® 80 Biomaster and Stomacher® 80 microBiomaster provide small tissue processing solutions for clinical, life science and microbiology applications. Capable of processing samples from 250µl to 80ml, the Biomaster range caters for the increased demand for smaller samples without reducing sample quality.

The Stomacher® 80 Biomaster (5ml to 80ml sample capacity), microBiomaster (250µl to 80ml sample capacity) and Stomacher® Bags allow samples to be processed safely and without risk of cross contamination while maximising visible cell recovery.

Tens of thousands of referenced methods are available and now with the enhanced sample capacity range of the microBiomaster the Stomacher® 80 is the World's most efficient and versatile sample preparation system for life science applications.

Choose the microBiomaster for increased sampling volume range with the new microBag for effective homogenization and sample recovery below 1ml volumes.

Why the Stomacher® 80?

Adjustable process times
& paddle speeds

Proven enhancement in
cell recovery & viability

No risk of cross
contamination

Process samples
from 250µl - 80ml

Cost effective sample
preparation




Stomacher® 80
Biomaster



Stomacher® 80
microBiomaster

Stomacher® 80 Biomaster



The Stomacher® **Bag** Range has been designed to fulfil all applications from blending, straining and storage including bags designed specifically for life science applications.

Enhanced Stomacher® **80** bag range features closure and strainer bag formats, double bags for biosafety and the new microBag for effective homogenisation and sample recovery below 1ml volumes.


Seward Bag racks provide the perfect storage for Stomacher® **Bag** sample handling, helping to provide a clean and productive working environment.

80 Bags

Product Code	Bag Type	Bag Size	Sample Size	Pack Quantity
BA6040	Standard Bag	105 x 155mm	5-80ml	2000 (sachets of 50)
BA6040CLR	Closure Bag	105 x 155mm	5-80ml	200 (sachets of 10)
BA6040/STR	Strainer Bag	105 x 155mm	5-80ml	200 (sachets of 10)
BA6040/STR/DBL	Strainer Double Bag	135 x 180mm	5-80ml	200 (sachets of 10)
BA6040/CLR/STR	Strainer Closure Bag	105 x 155mm	5-80ml	200 (sachets of 10)
BA6040/MICRO	Micro Bag	60 x 185mm	250µl-80ml	50

Accessories

Product Code	Description	Capacity	Dimensions	Weight
BA6098	Stomacher® 80 Bag Opener	1 Place	7 x 7 x 14cm	0.1kg
BA6090	Stomacher® 80 Bag Rack	10 Place	30 x 11 x 12cm	0.9kg
BA6051	Simplettes® Disposable straws	200 (8 sachets of 25 each)	24 x 18 x 18cm	3.7kg



www.seward.co.uk

Seward Limited
Dominion House
Easting Close
Worthing
West Sussex, BN14 8HQ
United Kingdom

Phone +44 (0)1903 823 077
Email info@seward.co.uk
Fax +44 (0)1903 219 233

Seward Laboratory Systems Inc. (USA)
574 NW Mercantile Place
Unit 107
Port Saint Lucie
FL 34986
USA

Phone + 772-621-8220
Fax + 772-621-8257
Email info@sewardusa.com

Seward International Pte Ltd
16 Raffles Quay
#33-03 Hong Leong Building
Singapore
048581

Phone +65 6519 2283
Email info@seward.co.uk