


Stomacher® 3500 Autoclavable blender bag

The first paddle blender bag designed to be autoclave resistant for disposal

This innovative technology aims to:

- Reduce single-use plastic consumption by eliminating the need to 'double-bag' the sample during autoclaving.
- Lower laboratory waste costs by having less materials to dispose of
- The same Seward Stomacher® 3500 autoclavable bag can be used to contain samples from blending to incubation to autoclaving for disposal



*Photo for illustrative purposes only

Good sample preparation is critical for dependable results

Seward makes products that offer real value to the laboratory. Seward Stomacher blenders have been proven over 40 years and are engineered to be reliable, strong, versatile, powerful and long-lasting. They are used to prepare over **8 million food samples in more than 130 countries** around the world and remain the best and most trusted piece of equipment in their field.



The gold standard in sample preparation

Home of the Stomacher®, Seward is the most highly referenced brand name in food safety and life science sample preparation in scientific publications. Focusing on transparency, accuracy and traceability, the Seward range helps to produce accurate, repeatable and robust results.

The Autoclavable Stomacher® bag

Following development and extensive evaluation by Seward, the new Stomacher® 3500 Autoclavable blender bag has been shown to be robust enough to withstand both the physical demands of the blending action and subsequent autoclaving. Compatible with the Seward Stomacher® 3500 the blender bags are suitable for environmental analysis, food microbiology, sterility testing and industrial mixing applications.

- ✓ Reduces plastic waste by 50% by eliminating the need to 'double-bag' samples for autoclaving.
- ✓ Made from high quality food grade material and extruded to 100 µm thickness, highly transparent and puncture resistant.
- ✓ Ideal for sample sizes ranging from 400-4500 ml.



Stomacher® 3500 Bags & Accessories

Product Code	Bag Type	Bag Size	Sample Size	Pack Quantity
BA6042	Standard Bag	380 x 510mm	400-3500ml	250 (sachets of 50)
BA6042/HD	Heavy Duty Bag	380 x 510mm	400-3500ml	200 (sachets of 50)
BA6042/ACL	Autoclavable Bag	390 x 520mm	400-4500ml	500 (sachets of 50)
BA6042/CLR	Closure Bag	380 x 510mm	400-4500ml	200 (sachets of 20)

Product Code	Description	Capacity	Dimensions	Weight
BA6094	Stomacher® 3500 Bag Rack	3 Place	42 x 15 x 35cm	5kg
BA6097	Stomacher® 3500 Bag Clips	N/A	N/A	N/A

Seward Limited
Technology Centre
Easting Close
Worthing
West Sussex, BN14 8HQ
United Kingdom

Seward Laboratory Systems Inc. (USA)
155 Keyland Court
Bohemia
NY 11716
USA

FermionX GmbH
Münstertäler Straße 12
Eschbach
DE-79427
Germany

Phone +44 (0)1903 524 600
Email info@seward.co.uk

Phone +1 631 337 1808
Email info@sewardusa.com

Email info@seward.co.uk

Seward is a part of the FermionX Group.