

Stomacher® Bags - Closure Bags

Product Codes by Sampling Volume

5-80ml: **BA6040/CLR, BA6040/STR/CLR**

80-400ml: **BA6041/CLR, BA6141/CLR**

400-4500ml: **BA6042/CLR**


TECHNICAL SPECIFICATIONS

- Stomacher® blender bag, with addition of bag closure device, consisting of wire reinforced tape
- Bag made from virgin grade, Polyolefin complex. Flexible and transparent
- Resistant to freezing and high temperatures (from -40°C to 80°C / -40°F to 176°F)
- No contact between the sample and the blender during homogenisation of the sample
- Suitable for short term storage or incubation of samples
- Compatible with any paddle lab blender
- Suitable for food contact. Not designed to be used for food packaging
- Irradiated sterile in sealed sachets, with batch traceable certificate
- Storage: in a temperate and dry place, away from humidity and light
- Available in 80ml, 400ml and 3500ml

Stomacher® 80 Biomaster Closure Bags


P/N. **BA6040/CLR**

- Max Blending Volume: 80ml
- Recommended blending Volume: 5-50ml
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)



P/N. **BA6040/STR/CLR**

- Max Blending Volume: 80ml
- Recommended blending Volume: 5-50ml
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 100x155ml (65µm thickness)
- Additional strainer material, consisting of microfiltered Polyolefin complex.



Seward Ltd

Technology Centre
Easting Close
Worthing, BN14 8HQ
United Kingdom
info@seward.co.uk

Seward Lab Systems Inc

155 Keyland Court
Bohemia
NY 11716
USA
info@sewardusa.com


FermionX GmbH

Münstertäler Str. 12,
Eschbach
DE-79427
Germany
info@seward.co.uk

Stomacher® 400 Closure Bags


P/N. BA6041/CLR

- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Straight bottom seal
- Pack size: 250 Bags (25 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)



P/N. BA6141/CLR


- Max Blending Volume: 400ml
- Recommended blending Volume: 100-250ml
- Round bottom seal
- Pack size: 250 Bags (25 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)



Stomacher® 3500/4500 Closure Bags

P/N. BA6042/CLR

- Max Blending Volume: 4500ml
- Recommended blending Volume: 3750ml
- Straight bottom seal
- Pack size: 200 Bags (20 x 10 sachets)
- Bag Dimensions: 177x305ml (65µm thickness)



Stomacher® is registered trademark of Seward Ltd.

SEWARD is a part of the FermionX Group. FermionX Ltd reserves the right to change and improve specifications of its products without notice. Please contact info@seward.co.uk or visit www.seward.co.uk for more information.

Seward Ltd

Technology Centre
Easting Close
Worthing, BN14 8HQ
United Kingdom
info@seward.co.uk

Seward Lab Systems Inc

155 Keyland Court
Bohemia
NY 11716
USA
info@sewardusa.com

FermionX GmbH

Münstertäler Str. 12,
Eschbach
DE-79427
Germany
info@seward.co.uk