



seward

World leaders in sample preparation

Introducing the new
range of Simplettes®

NEW ADDITION

BA6053 Simplette® Straws

Benefits

- Fits Socorex Acura® 810
- New clear packaging sachets with zip lock seal
- Cheap alternative to graduated pipette tips
- Straw length: 177mm, internal diameter 3.6-3.7mm
- Irradiated sterile

Other Simplette® Straws available
see over for more information

www.seward.co.uk

Other Simplette® straws available

Product Code	Description	Quantity	Internal Diameter
BA6051	Polypropylene Disposable Simplette® Straws -177mm	2000 per pack	4.1-4.2mm (Suitable for Barky Ultipette®)
BA6051/B	Blue Polypropylene Disposable Simplette® Straws - 177mm	2000 per pack	4.1-4.2mm (Suitable for Barky Ultipette®)
BA6052	Polypropylene Disposable Simplette® Straws - 250mm	1000 per pack	4.1-4.2mm (Suitable for Barky Ultipette®)
BA6053	Polypropylene Disposable Simplette® Straws – 177mm	2000 per pack	3.6-3.7mm (Suitable for Socorex Acura® 810)

For further information and details regarding Simplettes® or any Seward products please contact your local Seward distributor or visit: www.seward.co.uk

Seward Limited
Technology Centre
Easting Close
Worthing
West Sussex, BN14 8HQ
United Kingdom

Seward Laboratory Systems
Inc. (USA)
155 Keyland Court
Bohemia
NY 11716
USA

FermionX GmbH
Münstertäler Straße 12
Eschbach
DE-79427
Germany

Phone +44 (0)1903 524 600
Email info@seward.co.uk

Phone +1 631 337 1808
Email info@sewardusa.com Email info@seward.co.uk

Socorex Acura® 810 is a registered trademark of Socorex, Switzerland.
Barky Ultipette® is a registered trademark of Barky, England.