

Insterprep®

Ready-to-mix granulated media in a bag

Instantly soluble and sterile pre-measured media for fast sample preparation

A fast way to mix media ready for use instantly

Our Insterprep® bags and sachets allow laboratories to respond quickly to varying daily workflows, reducing waste, and saving preparation time.

Our Insterprep® range of pre-filled granular media bags and sachets contain exactly the right amount of dehydrated media for sample preparation.

Available in Buffered Peptone Water and Fraser Broth.

Just add sterile water and your sample.

The action of your Stomacher® paddle blender ensures the media is prepared at the same time as processing your sample, reducing a testing cycle from hours to minutes.

Choose from pre-filled and weighed Stomacher® sample bags or sachets that are ready to use with your own sample bags.

Insterprep is manufactured to ISO standards and is suitable for TVC and pathogen testing. Designed to work with our Stomacher® 400 Circulator, our Seward Stomacher® 400 EVO, Stomacher® Eco bags, and all other paddle type blenders.



Why Insterprep®?

Respond quickly to sample test requests

No more preparing of batches of new media at the end of the day/shift


Reduce waste by only preparing the media you need

Ease workflow bottlenecks by cutting media preparation time

Simple to prepare, no cleaning of media containers at the end of the day

	No Waste	Low odour/ powder residue	Easy to ship and store	Sample prepared in under 1 min
Powdered media			✓	
Prepared liquid media		✓		✓
Other granulated media sachets	✓	✓	✓	
Insterprep®	✓✓	✓✓	✓✓	✓✓*

*Insterprep media bags. Average time from opening bag, weighing sample and adding sterile water



Insterprep®

Ready-to-mix granulated media in a bag

- ✓ Flexibility to respond to varying workflow
- ✓ Easier to store compared with traditional bottles of media, saving storage space
- ✓ Instant preparation time
- ✓ Long shelf life for Insterprep, 2 years versus bulk prepared media bottles
- ✓ Reduce supply chain carbon emissions, shipping costs and waste disposal costs, compared to pre-mixed diluents
- ✓ Smooths laboratory workflow with no bottlenecks
- ✓ For all sizes of pre-enrichment

Save
space

Save
time

Reduce
waste

Insterprep®

Product		Product Code
Buffered Peptone Water (Bag, Granules for 90ml)	30 sachets of 10 bags	INS-2101-300NO
Buffered Peptone Water (Bag, Granules for 90ml)	6 Sachets of 10 bags	INS-2101-60NO
Buffered Peptone Water (Bag, Granules for 225ml)	30 Sachets of 10 bags	INS-2100-300NO
Buffered Peptone Water (Bag, Granules for 225ml)	6 Sachets of 10 bags	INS-2100-60NO
Buffered Peptone Water (Sachet, Granules for 90ml)	100 sachets	INS-2200-100NO
Buffered Peptone Water (Sachet, Granules for 3375ml)	250 sachets	INS-2300-250NO
Fraser Broth (Bag, Granules for 225ml)	30 sachets of 10 bags	INS-2102-300NO
Fraser Broth (Bag, Granules for 225ml)	8 Sachets of 10 bags	INS-2102-80NO
Fraser Broth (Sachet, Granules for 225ml)	100 sachets	INS-2202-100NO
Fraser Broth (Sachet, Granules for 1125)	50 sachets	INS-2302-50NO

Seward Limited
Technology Centre
Easting Close
Worthing
West Sussex, BN14 8HQ
United Kingdom

Phone +44 (0)1903 524 600
Email info@seward.co.uk

Seward Laboratory Systems Inc. (USA)
155 Keyland Court
Bohemia
NY 11716
USA

Phone +1 631 337 1808
Email info@sewardusa.com

FermionX GmbH
Münstertäler Straße 12
Eschbach
DE-79427
Germany

Email info@seward.co.uk