



# Stomacher® 4500 BAM



The robust Stomacher® 4500 is ideal for the BAM large volume testing procedures.

Expanding our range of blenders focused on food laboratories requirements, the Stomacher® 4500 BAM is the perfect tool for extraction from 375g composite samples.

The Stomacher® technique has become the world standard in sample preparation. The heavy duty Stomacher® 4500 BAM complements the Stomacher® 400 Circulator

and provides food laboratories with a high performance, high quality blender capable of homogenizing samples up to 4500ml.

### Why the Stomacher® 4500 BAM?

Ideal for 375g samples

Efficient Large Volume Stomacher®

Volume sample size:  
400 – 4500ml

Programmable Settings

3 year Warranty

For further information on Seward products please contact your local Seward distributor.

Seward Limited  
 Technology Centre  
 Easting Close  
 Worthing  
 West Sussex, BN14 8HQ  
 United Kingdom




Phone +44 (0)1903 524 600  
 Email info@seward.co.uk

Seward Laboratory Systems Inc. (USA)  
 155 Keyland Court  
 Bohemia  
 NY 11716  
 USA

Phone +1 631 337 1808  
 Email info@sewardusa.com

FermionX GmbH  
 Müntertäler Straße 12  
 Eschbach  
 DE-79427  
 Germany

Email info@seward.co.uk

Machine Model	Paddle Type	Sample weight*	Diluent Volume*
 4500/001/EU(240V) 4500/001/AJ/(110V)		<b>375g</b>	<b>3375ml</b>
Bags/Accessories	Standard Bags	Bag Rack	Bag Clip
	BA6042 BA6042/HD BA6042/ACL	BA6094	BA6097

\* as per the FDA BAM protocol protocol listed 20 September 2014 :  
[www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods](http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods)